

URIT

Advanced dual-magnetic circuit bead method
2-channel random & reliable hemostasis

URIT-600

Coagulation Analyzer



Global Diagnostics Supplier

URIT-600

Coagulation Analyzer

Specifications

Principle: dual-magnetic circuit bead method

Parameters:

PT (Prothrombin Time)

APTT (Activated Partial Thromboplastin Time)

TT (Thrombin Time) FIB (Fibrinogen)

Protein C, Protein S, HEP (High-molecular-weight Heparin),
LMWH(Low-molecular-weight Heparin),

LA(Lupus Anticoagulant),

Thrombin factor II, V, VII, X, VIII, IV, IX, X I, X II etc.

Detector: 2 channels

Incubator: 16 sample positions, 4 reagent positions (tunable incubation time)

Data storage: 500 test results

Display: 5.1" LCD screen

Print: Built in printer & external printer

Input / Output: RS-232 port and parallel interface

Ambient: work temperature: 15~35℃;
relative humidity: ≤80%RH

Power: AC100V~240V, 50/60Hz, 45W

Dimension: 380 mm(W)×340mm(D)×130mm(H)

Weight: 5.9 kg



Features

- 2 detectors up to 2 different parameters analysis
- Accurate incubator temperature
- Advanced method to eliminate interference from hemolysis, chyle, icterus, turbidity, plasma viscosity etc.
- Mini-reagent consumption, open reagent system
- Internal thermal printer with reference value range

URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Add: No. 4 EAST ALLEY, JIUHUA ROAD, GUILIN, GUANGXI 541001, P.R. CHINA

Tel: +86 773 2288583/2288589 Fax: +86 773 2288560/2288559

E-mail: export@uritest.com service@uritest.com Http://www.urit.com



Rev: 201010-4 (2010000)